



Deciding to select a better cleaning solution

Author

Ms. Nora Abdullah, CSSD Manager, Healthpoint Abu Dhabi, UAE
RN, Perioperative (MOH-Malaysia); Bachelor of Nursing (Monash-Australia);
CRCST, CIS, CHL (IAHCSMM-US)

Learning Objectives



1. Identify some of the critical selection criteria
2. Acquired, learned and understood the crucial factors
3. Share success and temporary defeats with the CS Communities
4. Profit by the experience
5. Invest concentration in improving the practices
6. Achieved the best sustainable outcomes

How to make better selection decisions?





1. A solution that is suitable for manual and mechanical cleaning



SUITABLE



2. Perform a comparative study between a few types of chemical




3. Washer-disinfector – two chemicals instead of three




4. A Detergent with disinfection



- a. Protect CS team from contagious materials
- b. For Devices incompatible with mechanical cleaning
- c. Specifically, tools exposed to infectious materials; MDRO, Hep. C
- d. Rapid contact time, broad-spectrum
- e. Allow immersion without damaging the tools



5. Acquiring history of the vendor consistent supplies to the end-users




6. Select the vendor that supply to more than three facilities



7. Vendor after sales support



- Assist continuous monitoring, educating the users, acquiring feedback, reporting



8. Establish a win-win relationship

- The Vendor/CSSD calculate
 - Suitable dose
 - Achieved hygienic condition
- Measure volume used per week/month
- Negotiate a term/agreement to supply for a specific duration



9. Establish a win-win relationship

- Vendor moves the reserved stocks
 - support increase & decrease demand
 - between facilities
- Create a relationship of trust and harmony



10. The return of investment-1



- Dosing pump calibrated
 - Manual, soaking, ultrasonic washing
- Complimentary dosing machine
- Regular inspection
 - supplier & CSSD
 - dosing pump, check active implementation
- Vendor assists monitoring expiry dates, over/under stock, no wastage



11. The return of investment-2



- Continuous evaluation
 - Optimum outcome
 - Opportunity to reduce the chemical concentration
- Monitoring
 - Water Quality, Chemical Efficacy
 - Apply minimum dose, Optimum hygiene
- Using two chemicals instead of three types; washer-disinfector



12. Conforming to the local regulatory agencies & protocols

- Product registration certificate with local health authority
- MOH, FDA, OSHA
- Instruction for Used
- Material Safety Data Sheet
- Certificate of material compatibilities to washer disinfectant, ultrasonic machine & device types.



13. References



- Seavey, R. Association for the Advancement of Medical Instrumentation. Sterile ... ANSI/AAMI Standards ST79:2017 Comprehensive guide to steam sterilization and sterility assurance in ... ST79:2017 Updates – Cleaning agents When using an automated chemical delivery, p.44.
- Central Service Technical Manual (CRCST) Workbook 8th Edition In-text: (Central Service Technical Manual (CRCST) Workbook 8th Edition, 2018) p.137.
- Reprocessing of Instruments to Retain Value In-text: (A-k-i.org, 2019) p.33. Available at: http://www.a-k-i.org/fileadmin/downloads/broschueren/rot/rb_gb_11.web.pdf [Accessed 10 Sep. 2019].

Healthpoint: Central Sterile Supply Department

Evaluation Checklist for RIMD (Reusable Invasive Medical Device) Cleaning Solution	
Select <input checked="" type="checkbox"/> the appropriate remarks and submit the documents to support this evaluation.	
1. Product Name.	a. b. c. d.
2. Supplier and Manufacturer Name.	
3. List the # and the Name of the Facility you are supplying in the country.	a. # of Facility: b. Facility Name:
4. The Chemical and dosage for a. Washer-Disinfectant	a. Chemical Name: Dosage for 27 liters of water per cycle: <input type="text"/> ML b. Chemical Name: Dosage for 27 liters of water per cycle: <input type="text"/> ML c. Chemical Name: Dosage for 27 liters of water per cycle: <input type="text"/> ML
b. Ultrasonic Bath	Chemical Name: Dosage per MI/liter of water:
c. Manual Cleaning	Chemical Name: Chemical Dosage per MI/liter of water:
5. The Chemical added properties.	Chemical Name: <input type="checkbox"/> Disinfectant <input type="checkbox"/> Proteases <input type="checkbox"/> Amylases <input type="checkbox"/> Lipases <input type="checkbox"/> NA
6. The Chemical Compatibilites to materials and methods.	Chemical Name: <input type="checkbox"/> Washer-Disinfectant <input type="checkbox"/> Ultrasonic Bath <input type="checkbox"/> Manual Cleaning <input type="checkbox"/> Flexible Endoscope <input type="checkbox"/> Rigid Endoscope
7. Select the Historical undesirable events between Supplier and End-User.	<input type="checkbox"/> Shipment Delay <input type="checkbox"/> Manufacturer Delay <input type="checkbox"/> Supplier Out of Stock <input type="checkbox"/> User Over-stock & expired <input type="checkbox"/> Short Expiry Date <input type="checkbox"/> No Offer for Contingency <input type="checkbox"/> NA <input type="checkbox"/> Other _____ # of undesirable events: <input type="checkbox"/> 0 <input type="checkbox"/> 1 to 3 <input type="checkbox"/> 4 & above Between the period of: <input type="checkbox"/> 6 months <input type="checkbox"/> 1 year <input type="checkbox"/> 2 years <input type="checkbox"/> 3 years
8. Total cleaning cycle time to achieve hygenic condition and safe for handling.	a. Washer Disinfectant: <input type="text"/> minutes b. Ultrasonic Bath: <input type="text"/> minutes c. Manual Contact Time: <input type="text"/> minutes
9. Select one or more of the Product Label Statements.	<input type="checkbox"/> Rust Inhibiting <input type="checkbox"/> Odor Suppressant <input type="checkbox"/> Chloride Free <input type="checkbox"/> Disinfectant <input type="checkbox"/> Low Foam <input type="checkbox"/> Non-Corrosive <input type="checkbox"/> No Residual Stain <input type="checkbox"/> Compatible to Delicate RIMD (Ophthalmology, Micro devices) <input type="checkbox"/> Disposed in Common Sewage

14. The End

